




# GRASON & ASSOCIATES, LLC



P.O. Box 289, Rindge, NH 03461 USA

Tel. (603) 899-3089 Fax. (603) 899-9747 [www.grasonassociates.com](http://www.grasonassociates.com)

## What eartips fit my instrument?



Please use this guide to find out which eartips fit your instrument.

The list is sorted by instrument make. Some instruments may be listed twice due to differences in probe or other optional features. Be aware that some instruments or makes may have a note associated with them.

### To Search this Document:

Use the find function located under the edit menu or press ctrl-f.

It is recommended to search by manufacturer and look at multiple results, not just the first. Some instruments fit multiple types of eartips depending on probe.



If unsure of which eartip fits, please ask for a sample to try.

***Know of Other Instruments that single use eartips fit?***

Please let us know at: [sales@grasonassociates.com](mailto:sales@grasonassociates.com)

Document Format will be as below:

### Instrument Make

Instrument Name X ..... KR, RS, MF

Instrument Name Y ..... IA

## American Electromedics

Note: To our knowledge the KR, RS and MF series of eartips will fit the older American Electromedics Instruments, although fit is not exact. To Use KR, RS or MF eartips with an American Electromedics instrument, the eartips must be pushed onto the probe a long way until the forward edge of the probe is flush with the front of the eartip. The only known option for exact fit on these instruments is to find Maico Rock eartips.

81 .....	KR, RS, MF	85 .....	KR, RS, MF
83 .....	KR, RS, MF	86 .....	KR, RS, MF

## Amplaid

Note: To our knowledge the IT series of eartips will fit Amplaid Instruments.

## Amplivox

Otowave models with <b>White Probe</b> .....	KR, RS, MF	Otowave models with <b>Clear Probe</b> .....	IA
---	------------	---	----

## Danplex

Note: To our knowledge the KR, RS and MF series of eartips will fit Danplex Instruments.

## EAR

Note: To our knowledge the KR, RS and MF series of eartips will fit the EAR 3A Earphones with EAR3-06 adapter.

## Echodia

Elios .....	KR, RS, MF	Babyscreen .....	KR, RS, MF
Otodia .....	KR, RS, MF		

## Etymotic Research

Ero Scan DPOAE .....	KR, RS, MF	ER10D .....	KR, RS, MF
ER10 Low Noise™		ER3A Earphone with probe	
Microphone ER10B .....	KR, RS, MF	tip adapter .....	KR, RS, MF

## Grason Stadler

Tympstar .....	KR, RS, MF	GSI 20 through 39 .....	KR, RS, MF
Tympstar Pro .....	IA	1720 through 1739 .....	KR, RS, MF
Audera .....	GS	Corti .....	IA
Audio Screener .....	GS	70 .....	GS

## Intelligent Hearing Systems

Note: To our knowledge the KR, RS and MF series of eartips will fit the Smart OAE with 10D probe.

---

### Interacoustics

Note: To our knowledge many Interacoustics instruments can fit two different eartips depending on the model vintage. As a general rule for AT235, AT235h, AA220, AA222, and MT10: clear probes use the IA eartip, while white, grey or black probes use the KR, RS, MF eartips. If unsure, please ask us or your distributor for a sample to check fit.

Otoread .....	KR, RS, MF	Eclipse OAE (NEOPT probe) .....	IA
Titan .....	IA	Eclipse (clear probe) .....	IA
AT235 (clear probe) .....	IA	IMP440 .....	IA
AT235 .....	KR, RS, MF	MT10 (clear probe) .....	IA
AA22X (clear probe) .....	IA	MT10 .....	KR, RS, MF
AA22X .....	KR, RS, MF	MT5 .....	IA
BET50 compatible .....	KR, RS, MF	MTP .....	IA
BET55 compatible .....	IA	EP25 .....	IA

---

### Inventis SRL

Flute .....	IT	Viola .....	IT
Clarinet .....	IT		

---

### Kamplex

Note: To our knowledge the KR, RS and MF series of eartips will fit Kamplex Instruments.

---

### Labat Srl

Echolab .....	KR, RS, MF		
---------------	------------	--	--

---

### Madsen

Otoflex .....	MO	Alpha .....	PM
Zodiac .....	IA	Capella 2 .....	KR, RS, MF
Zodiac 901 .....	KR, RS, MF	Accuscreen .....	IA

---

## Maico

Note: To our knowledge many Maico instruments can fit two different eartips depending on the model vintage. As a general rule for Otowave, EasyTymp, MI22, MI23, MI24, MI26, MI34, MI44: clear probes use the IA eartip, while white, grey or black probes use the KR, RS, MF eartips. If unsure, please ask us or your distributor for a sample to check fit.

Otowave (clear probe) .....	IA	EasyTymp (clear probe) .....	IA
Otowave .....	KR, RS, MF	EasyTymp .....	KR, RS, MF
MA790 .....	KR, RS, MF	630 .....	KR, RS, MF
MA800 .....	KR, RS, MF	MI 22 (clear probe) .....	IA
Eroscan .....	KR, RS, MF	MI 22 .....	KR, RS, MF
Eroscan Pro .....	KR, RS, MF	Race Car Tymp (clear probe) .....	IA
Eroscan Pro + Tymp .....	KR, RS, MF	Race Car Tymp .....	KR, RS, MF
Eroscan Plus DPOAE .....	IA	MI 24 (clear probe) .....	IA
MI 23 (clear probe) .....	IA	MI 24 .....	KR, RS, MF
MI 23 .....	KR, RS, MF	MI 34 (clear probe) .....	IA
MI 26 (clear probe) .....	IA	MI 34 .....	KR, RS, MF
MI 26 .....	KR, RS, MF	MI 44 (clear probe) .....	IA
620 .....	KR, RS, MF	MI 44 .....	KR, RS, MF

## MedRx

Otowave (clear probe) ..... IA

## Micro Audiometrics

Earscan Acoustic Impedance Audiometer ..... KR, RS, MF

## Neurosoft

Neuro-Audio-Screen ..... KR, RS, MF | Neuro-Audio/OAE ..... KR, RS, MF

## Oscilla

TSM 300 ..... KR, RS, MF | TSM 400 ..... KR, RS, MF

## Otodynamics

Note: Grason & Associates sells two types of eartips that to our knowledge fit the Otodynamics instruments. TEOAE only probes that look like ☺ on the end will fit our OT0XX and OT2XX eartips. TEOAE + DPOAE probes that look like ☹ on the end will fit our OT3XX eartips.

ECP probe .....	OT0XX, OT2XX	UGD probe .....	OT3XX
UGS probe .....	OT0XX, OT2XX	SND probe .....	OT3XX
SNS probe .....	OT0XX, OT2XX	SGD probe .....	OT3XX
SGS probe .....	OT0XX, OT2XX		

## Path Medical

Sentiero .....	PM
----------------	----

## Resonance

Note: Eartip fit for the model numbers listed below will apply to all optional models. As an example, a listing for r25c will also fit the r25c-pt.

r16m .....	IT	r15c .....	IT
r26m .....	IT	r25c .....	IT
r36m .....	IT	r35c .....	IT

## Siemens

Note: To our knowledge the KR, RS and MF series of eartips will fit Siemens Instruments.

## Virtual

Note: To our knowledge the KR, RS and MF series of eartips will fit Virtual Instruments.

## Vivosonic

Integrity OAE .....	KR, RS, MF
---------------------	------------

## Welch Allyn

262 .....	KR, RS, MF	Micro Tymp .....	WA
39500 OAE Screener .....	KR, RS, MF	Micro Tymp 2 .....	WA