

POLY- SPECTRUM- 8/V

digital ECG system
for small and medium-sized
animals




FEATURES

- Recording ECG from standard limb leads (I, II, III), augmented limb leads (aVL, aVR, aVF) and precordial leads (V1...V6).
- Skin-safe electrodes that prevent the risk of infection due to skin injury and ensure better skin contact if compared to the needle electrodes.
- Automatic QRS detection and calculation of all temporary parameters using animal-specific algorithm (for the dog, cat, etc.)
- Automatic ECG measurement and interpretation, and report generation.
- Can be used as ECG monitor during surgery.

EXAMINATIONS




ECG of the dog



ECG of the dog. Averaged QRS complexes

DELIVERY SET



- Poly-Spectrum-8/V electronic unit
- 4-lead ECG cable (2 m)
- Veterinary ECG electrode with alligator clip (8 pcs.)
- Electrode gel Unimax (250 g)
- Technical and user manuals
- Software
- Transportation bag